

Data Collection & Data Analysis

Basic Processing with Mosflm


Gwyndaf Evans

Diamond Light Source, UK.

Contents

- some thoughts on data collection
 - simple processing with *Mosflm*
 - data collection strategy with *Mosflm*
-

MOSFLM and SCALA in CCP4



Optimization of Data Collection

- Pre-process at least one image (preferably two at 90° to each other) to obtain:
 - Cell parameters, crystal orientation and putative Laue group
 - Estimate of mosaicity
 - Effective resolution limit
 - Crystal to detector distance
 - Exposure time
 - Strategy for data collection
 - **Remember!**
 - This is the last experimental stage - if you collect bad data now you are stuck with it.
 - No data processing program can rescue the irredeemable!
-

What is needed prior to running Mosflm?

- ❑ X-ray images
 - ❑ Experimental details (e.g. detector type, direct beam position, wavelength, etc)
 - ❑ The program itself and a computer to run it on!
-

[localhost:~/test/muldlx1] harry% ipmosflm

***** Version 6.2.3 for Image plate and CCD data 10th July 2003 *****

A. G. W. Leslie, MRC Laboratory Of Molecular Biology, HILLS ROAD, CAMBRIDGE CB2 2QH, UK

E-mail andrew@mrc-lmb.cam.ac.uk

New auto-indexing using DPS due to Ingo Steller Robert Bolotovskiy and Michael Rossmann

(1998) J. Appl. Cryst. 30, 1036-1040

Original auto-indexing using REFIX due to Wolfgang Kabsch (Kabsch, W. (1993),

J. Appl. Cryst. 24, 795-800.)

X-windows interface using xdl_view due to John Campbell (Daresbury Laboratory, UK.)

(Campbell, J. W. (1995) J. Appl. Cryst. 28, 236-242.

MOSFLM => image muldlx1_301.mar2000

MOSFLM => go

(Q) QOPEN: file opened on unit 1 Status: READONLY

<!--SUMMARY_BEGIN-->

Logical Name: muldlx1_301.mar2000 Filename: muldlx1_301.mar2000

<!--SUMMARY_END-->

Crystal to detector distance of 250.00 mm taken from image header

Wavelength of 1.54180A taken from image header

Crystal to detector distance of 250.00 mm taken from image header

Wavelength of 1.54180Å taken from image header

Pixel size of 0.1500 mm taken from image header.

Start and end phi values for image 1 from image header are 279.00 and 280.00 degrees.
image FILENAME: muldxx1_301.mar2000

```
(Q) QOPEN: file opened on unit 1    Status: READONLY
<B><FONT COLOR="#FF0000"><!--SUMMARY_BEGIN-->
  Logical Name: muldxx1_301.mar2000    Filename: muldxx1_301.mar2000
<!--SUMMARY_END--></FONT></B>
```

The red circle denotes the region behind the backstop shadow
(Use BACKSTOP keyword to set this.)

Edits allowed

Processing params

```

a      : 0.00
b      : 0.00
c      : 0.00
alpha  : 0.00
beta   : 0.00
gamma  : 0.00
PsiX   : 0.00
PsiY   : 0.00
PsiZ   : 0.00
Mosaic : 0.000
Divh   : 0.000
Divv   : 0.000
Lambda : 1.542
Distance: 250.00
Beam X : 149.60
       Y : 149.80
CCOMEGA : 0.000
ROFF   : 0.00
TOFF   : 0.00
YSCAL  : 1.0000
Pick area: X: 11
         Y: 11
Int threshold: 20
Vector scale 1
Two theta 0.00
Resolution 0.00
*SPOT SEARCH*
Threshold 5.00
Rmin 15.00
Rmax 135.00
X offset 0.00
Y offset 0.00
Min X size 0.50
Max X size 2.00
Min Y size 0.50
Max Y size 2.00
Min no of pix 6
X splitting 0.30
Y splitting 0.30
*AUTOINDEXING*
Min I/sig(I): 20

```

Prompts On
 Update display: No
 After refinement No
 After integration No
 Timeout mode Off

 Select item

Main menu

```

Read image
Find spots
Edit spots
Clear spots
Select images
Autoindex
Estimate mosaicity
Predict
Clear prediction
Adjust
Refine cell
Integrate
Strategy
Keyword input
Find hkl
Pick
Measure cell
Circles
Beam / backstop
Save/Exit

```

Output

```


Pixel X,Y 2001 0
XC,YC mm 300.0 0.0
Resolution 0.00
Indices 0 0 0
F Phi 0.00 width 0.00
Intensity 0
Sigma 0

Spacing A 0.000
Average 0.0
Rms 0.0
Number 0
Zoomfactor 0
Circle resolution A
0.0 0.0 0.0 0.0
Phi 279.00 280.00
Missets ThetaX,Y,Z
0.00 0.00 0.00

```

 muldx1_301.mar2000

Min 1 Max 2207 Cursor position
 Overlay on Contrast
 Colour Black on whi Mag x4 PS Zoom



Blue: fulls, Yellow: partials, Red: overlaps
 Green: too wide in phi

Processing params

```

a      :      0.00
b      :      0.00
c      :      0.00
alpha  :      0.00
beta   :
gamma  :

```

Main menu

Read image

Find spots

Min 1 Max 2207 Cursor position

Overlay onContrast Colour Black on whi Mag x4

xdl_io_window

 Input reply

```

Find spots for autoindexing. Parameters determining
the spot finding are listed under *SPOT SEARCH*
in the "processing params" menu. Pixels with values
greater than "Threshold" sds above background are
considered to belong to spots. The program will set
a suitable threshold automatically, but for
"diffficult" images it can be set manually.
The spot search will be between Rmin and Rmax mm of
the direct beam position.
The median spot size (in X and Y) will be determined
and Xmin,Xmax,Ymin,Ymax are the limits on the spot
sizes as a function of this median size. Reduce
the max values if spots are not well resolved.
"Min no of pix" is the minimum number of pixels in
in a spot. Use the "Min I/sig(I)" parameter to test
the effect of changing the intensity threshold
applied in the actual autoindexing.

```

Do you wish to continue ? (Y):

Do you want to find spots manually ? (N) _

```

Int
Vect
Two
Reso
*SPD
Thre
Rmin
Rmax
X of
Y of
Min
Max
Min
Max
Min
X sp
Y sp
*AUT
Min


```

```

Prom
Upda
Afte
Afte

```

Time



```

Missets ThetaX,Y,Z
0.00 0.00 0.00

```

Edits allowed

Processing params

```

a      : 126.19
b      : 126.19
c      : 74.31
alpha  : 90.00
beta   : 90.00
gamma  : 120.00
PsiX   : 0.00
PsiY   : 0.00
PsiZ   : 0.00
Mosaic : 0.510
Divh   : 0.000
Divv   : 0.000
Lambda : 1.542
Distance: 250.00
Beam X : 149.69
      Y : 149.71
CCOMEGA : 0.000
ROFF   : 0.00
TOFF   : 0.00
YSCAL  : 1.0000
Pick area: X: 11
      Y: 11
Int threshold: 20
Vector scale 1
Two theta 0.00
Resolution 2.88
*SPOT SEARCH*
Threshold 4.92
Rmin 15.00
Rmax 135.00
X offset 0.00
Y offset 0.00
Min X size 0.50
Max X size 2.00
Min Y size 0.50
Max Y size 2.00
Min no of pix 6
X splitting 0.30
Y splitting 0.30
*AUTOINDEXING*
Min I/sig(I): 20

Prompts 
Update display: 
After refinement 
After integration 

Timeout mode 
    
```

Select item

Main menu

- Read image
- Find spots
- Edit spots
- Clear spots
- Select images
- Autoindex
- Estimate mosaicity
- Predict
- Clear prediction
- Adjust
- Refine cell
- Integrate
- Strategy
- Keyword input
- Find hkl
- Pick
- Measure cell
- Circles
- Beam / backstop
- Save/Exit

Output

```

Pixel X,Y 2001 0
XC,YC mm 300.0 0.0
Resolution 0.00
Indices 0 0 0
F Phi 0.00 width 0.00
Intensity 0
Sigma 0


Spacing A 0.000
Average 0.0
Rms 0.0
Number 0
Zoomfactor 0
Circle resolution A
0.0 0.0 0.0 0.0
Phi 279.00 280.00
Missets ThetaX,Y,Z
0.00 0.00 0.00
    
```

muld1x1_301.mar2000

Min 1 Max 2207 Cursor position

Overlay on Contrast

Colour Black on whi Mag x4



Blue: fulls, Yellow: partials, Red: overlaps
Green: too wide in phi

xdl_io_window

Input reply

(see terminal window for full output)
Autoindexing will use spots from image(s): 301

The new auto-indexing option works well with either single or multiple images, and with 100 - 750 spots: it is particularly good if your cell is very anisotropic or if you have strong diffraction.

N.B. it tries to estimate the longest cell edge. but you can override this if you wish.

Do you want to try the new auto-indexing? (Y):

Do you want to fix the detector distance (Y) ? :
Filename for final orientation matrix (muld1x1_301.mat):

Maximum expected cell edge (Angstroms) [162]:

Do you want to pre-refine the solutions? (N):

Do you want to proceed (Y):_

tion

Zoom

xdl_io_window

Input reply

Only solutions with PENALTY less than 200 are listed, a complete list is given in the terminal window

No	PENALTY	LATT	a	b	c	alpha	beta	gamma	Possible spacegroups
15	200	mP	74.34	126.23	125.97	119.9	90.2	90.0	P2, P21
14	139	mC	262.84	74.34	126.23	90.0	118.8	73.8	C2
13	139	oI	74.34	125.97	230.86	90.1	108.6	89.8	I222, I212121,
12	138	mC	262.40	74.34	126.23	90.0	118.6	73.8	C2
11	138	mI	74.34	230.86	125.97	90.1	90.2	71.4	C2
10	3	hP	125.97	126.23	74.34	90.0	90.2	119.9	P3, P31, P32, P312, P321, P3112, P3121, P3212, P3221 P6, P61, P65, P62, P64, P63, P622, P6122, P6522, P6222, P6422, P6322 C222, C2221
9	3	oC	126.22	218.35	74.34	89.9	90.3	89.9	C2
8	3	mC	126.22	218.35	74.34	89.9	90.3	89.9	C2
7	2	mC	218.32	126.23	74.34	90.0	90.3	90.1	C2
6	1	mC	125.97	218.77	74.34	89.8	90.2	90.0	C2
5	1	mC	218.77	125.97	74.34	90.2	90.2	90.0	C2
4	1	mP	125.97	74.34	126.23	90.0	119.9	90.2	P2, P21
3	1	oC	125.97	218.77	74.34	89.8	90.2	90.0	C222, C2221
2	0	aP	74.34	125.97	126.22	60.1	89.7	89.8	P1
1	0	aP	74.34	125.97	126.23	119.9	90.0	90.2	P1

Select a solution AND a spacegroup from list above (eg 3 p42) or 0 to abandon or T to change min I/sig(I):10 p64

The solution and direct beam position will now be refined; reflections which deviate by more than the sigma cutoff from their calculated position will be excluded from the refinement.

Positional sigma cutoff [2.50]:
Refining solution #10 with P64 (number 172) symmetry imposed

Using 146 indexed reflections (out of 155 spots found, ($\Delta(XY) \leq 2.5 \text{ sigma}$)),
final sd in spot positions is 0.12mm and in phi 0.23 degrees
Refined cell parameters 126.19 126.19 74.31 90.00 90.00 120.00

Do you want to update cell parameters (Y):


Beam coordinates of 149.60 149.80 have been refined to 149.69 149.71

This is a shift of 0.12mm or 0.035 times the minimum spot separation of ca 3.53mm.
Do you want to accept the new direct beam position? (Y) :

Do you want to accept this solution (Y) :_

Processing params		Main menu	
a	: 0.00	Read image	
b	: 0.00	Find spots	
c	: 0.00	Edit spots	
alpha	: 0.00	Clear spots	
beta	: 0.00	Select images	
gamma	: 0.00	Autoindex	
PsiX	: 0.00	Estimate mosaicity	
PsiY	: 0.00	Predict	
PsiZ	: 0.00	Clear prediction	
Mosaic	: 0.000	Adjust	
Divh	: 0.000		
Divv	: 0.000		
Lambda	: 1.542		
Distance:	250.00		
Beam X	: 149.60		
Y	: 149.80		

Min 1 Max 2207 Cursor position
 Overlay on Contrast
 Colour Black on whi Mag x4



xdl io window

 Input reply

Only solutions with PENALTY less than 200 are listed, a complete list is given in the terminal window

No	PENALTY	LATT	a	b	c	alpha	beta	gamma	Possible spacegroups
15	200	mP	74.34	126.23	125.97	119.9	90.2	90.0	P2, P21
14	139	mC	262.84	74.34	126.23	90.0	118.8	73.8	C2
13	139	oI	74.34	125.97	230.86	90.1	108.6	89.8	I222, I212121,
12	138	mC	262.40	74.34	126.23	90.0	118.6	73.8	C2
11	138	mI	74.34	230.86	125.97	90.1	90.2	71.4	C2
10	3	hP	125.97	126.23	74.34	90.0	90.2	119.9	P3, P31, P32, P312, P321, P3112, P3121, P3212, P3221 P6, P61, P65, P62, P64, P63, P622, P6122, P6522, P6222, P6422, P6322 C22, C2221
9	3	oC	126.22	218.35	74.34	89.9	90.3	89.9	C2
8	3	mC	126.22	218.35	74.34	89.9	90.3	89.9	C2
7	2	mC	218.32	126.23	74.34	90.0	90.3	90.1	C2
6	1	mC	125.97	218.77	74.34	89.8	90.2	90.0	C2
5	1	mC	218.77	125.97	74.34	90.2	90.2	90.0	C2
4	1	mP	125.97	74.34	126.23	90.0	119.9	90.2	P2, P21
3	1	oC	125.97	218.77	74.34	89.8	90.2	90.0	C222, C2221
2	0	aP	74.34	125.97	126.22	60.1	89.7	89.8	P1
1	0	aP	74.34	125.97	126.23	119.9	90.0	90.2	P1

Select a solution AND a spacegroup from list above (eg 3 p42) or 0 to abandon or T to change min I/sig(I):10 p64

The solution and direct beam position will now be refined; reflections which deviate by more than the sigma cutoff from their calculated position will be excluded from the refinement.

Positional sigma cutoff [2.50]:

Refining solution #10 with P64 (number 172) symmetry imposed

Using 146 indexed reflections (out of 155 spots found, {delta(XY) <= 2.5 sigma}),

final sd in spot positions is 0.12mm and in phi 0.23 degrees

Refined cell parameters 126.19 126.19 74.31 90.00 90.00 120.00

Do you want to update cell parameters (Y):

Beam coordinates of 149.60 149.80 have been refined to 149.69 149.71

This is a shift of 0.12mm or 0.035 times the minimum spot separation of ca 3.53mm.

Do you want to accept the new direct beam position? (Y) :

Do you want to accept this solution (Y) :_

Edits allowed

Processing params

```

a      : 126.19
b      : 126.19
c      : 74.31
alpha  : 90.00
beta   :
gamma  :
PsiX   :
PsiY   :
PsiZ   :
Mosaic :
Divh   :
Divv   :
Lambda :
Distance:
Beam X :
      Y :
CCOMEGA :
ROFF   :
TOFF   :
YSCAL  :
Pick area: X:
      Y:
Int threshold:
Vector scale
Two theta
Resolution
*SPOT SEARCH*
Threshold
Rmin
Rmax
X offset
Y offset
Min X size
Max X size
Min Y size
Max Y size
Min no of pix
X splitting
Y splitting
*AUTOINDEXING*
Min I/sig(I):

Prompts
Update display:
After refinement
After integration

Timeout mode Off
    
```

Select item

Main menu

Read image

muldlx1_301.mar2000

Min 1 Max 2207 Cursor position

Overlay on Contrast

Colour Black on whi Mag x4 **PS** **Zoom**

xdl_io_window

Input reply

For this spacegroup one segment is usually sufficient but two segments may give improved accuracy.

Give number of segments (1) :

Image number for first image of segment 1 (301) :

Image identifier (muldlx1) :

Use phi values from image header ? (Y):

Number of images in this segment (2) :

Use the current crystal orientation (Y) :


Filename for final orientation matrix (muldlx1_301.mat):

Post refining cell using 1 segments

Segment 1 images 301 to 302

Image identifier muldlx1

Do you want to proceed (Y):_



```

Number      0
Zoomfactor  0
Circle resolution A
      0.0  0.0  0.0  0.0
Phi      279.00  280.00
Missets ThetaX,Y,Z
      0.00  0.00  0.00
    
```

Blue: fulls, Yellow: partials, Red: overlaps
Green: too wide in phi

MOSFLM Image Display

Edits allowed

Processing params

```

a      : 126.19
b      : 126.19
c      : 74.31
alpha  : 90.00
beta   : 90.00
gamma  : 120.00
PsiX   : 0.00
PsiY   : 0.00
PsiZ   : 0.00
Mosaic : 0.510
Divh   : 0.000
Divv   : 0.000
Lambda : 1.542
Distance: 250.00
Beam X : 149.69
      Y : 149.71
CCOMEGA : 0.000
ROFF   : 0.00
TOFF   : 0.00
YSCAL  : 1.0000
Pick area: X: 11
      Y: 11
Int threshold: 20
Vector scale 1
Two theta 0.00
Resolution 2.88
*SPOT SEARCH*
Threshold 4.92
Rmin 15.00
Rmax 135.00
X offset 0.00
Y offset 0.00
Min X size 0.50
Max X size 2.00
Min Y size 0.50
Max Y size 2.00
Min no of pix 6
X splitting 0.30
Y splitting 0.30
*AUTOINDEXING*
Min I/sig(I): 20

Prompts  On
Update display:  No
After refinement  No
After integration  No

Timeout mode  Off
    
```

Select item

Main menu

-
-
-
-
-
-
-
-
-
-
-
-
-

Output

```

Pixel X,Y 2001 0
XC,YC mm 300.0 0.0
Resolution 0.00
Indices 0 0 0
F Phi 0.00 width 0.00
Intensity 0
Sigma 0


Spacing A 0.000
Average 0.0
Rms 0.0
Number 0
Zoomfactor 0
Circle resolution A
0.0 0.0 0.0 0.0
Phi 279.00 280.00
Missets ThetaX,Y,Z
0.00 0.00 0.00
    
```

muldlx1_301.mar2000

Min 1 Max 2207 Cursor position

Overlay on Contrast

Colour Black on whi Mag x4



xdl_io_window

Central profile:muldlx1_301.mar2000

```

-1-1-1-1-1-0-1-1-0-1-1-1-1-0-1-1-1-1-0-0-1-1-1
-1-1-1-1-0-1-1-1-1-1-1-1-1-1-1-1-1-1-0-1-1-1-0
-1-1-1-1-1-1-1-1-1-1-1-0-1-1-1-1-1-1-1-1-0-0-1-1
-1-0-1-1-1-1-1-1-1-1-1-1-1-1-1-2-1-1-1-1-1-1-0
-0-1-0-0-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-0-1-1
-1-0-1-1-1-1-1-1-1 2 2 3 3 3 2 2*2*2-1-1-1-1-1-1
-1-1-1-1-1-1-1-1 2 2 3 5 6 6 5 4 3*2-1-1-1-1-1-1
-1-1-1-1-0-1 1 2 4 8 C E E B 7 5 3*2-1-1-1-1-1-1
-1-1-1-1-1-1 2 4 8 E M P N I C 7 3*2-1-1-1-1-1-0
-1-1-1-1-1-2 2 6 C M U X U N F 9 4*2*2-1-1-0-1-1
-1-1-1-1-1-1 3 8 G Q Y Z V O G 9 5*2-1-1-2-1-1-1
-1-1-1-1-0-1 3 8 G P V X S M F 8 4*2-2-1-1-1-1-1
-1-1-1-1-1-2 3 6 B I O P L H B 6 3*2-1-1-1-1-1-1
-1-1-1-1-1-1 2 3 6 A D E C A 6 4 2-1-1-1-0-1-1-1
-1-1-1-1-1-1 2 3 5 6 7 6 5 3 2-1-1-1-1-1-1-1-1
-1-1-1-1-1-1-1-1 2 3 3 3 3 3 2-1-1-1-1-1-1-1-0
-1-1-1-1-1-1-1-1-1-2-2*2-1-2-1-1-1-1-1-1-1-1-1-1
-0-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-0-1-1-1-0
-1-1-1-1-1-1-1-1-1-0-1-1-1-1-1-1-1-1-0-1-1-1-1
-1-0-1-1-0-0-1-1-1-1-1-1-1-1-1-0-0-1-1-1-1-1-1
-0-1-0-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-0-0-1-1-0
    
```

Rms residual 0.06mm, weighted residual 0.8

Blue: fulls, Yellow: partials, Red: overlaps
Green: too wide in phi

Edits allowed

Processing params

```

a      : 125.79
b      : 125.79
c      : 74.12
alpha  : 90.00
beta   : 90.00
gamma  : 120.00
PsiX   : 0.10
PsiY   : -0.01
PsiZ   : -0.01
Mosaic : 0.753
Divh   : 0.000
Divv   : 0.000
Lambda : 1.542
Distance:
Beam X :
      Y :
CCOMEGA :
ROFF   :
TOFF   :
YSCAL  :
Pick area: X:
          y:
Int threshold:
Vector scale
Two theta
Resolution
*SPOT SEARCH*
Threshold 4.91
Rmin      15.00
Rmax     135.00
X offset  0.00
Y offset  0.00
Min X size 0.50
Max X size 2.00
Min Y size 0.50
Max Y size 2.00
Min no of pix 6
X splitting 0.30
Y splitting 0.30
*AUTOINDEXING*
Min I/sig(I) 20

Prompts       On
Update display:
After refinement  No
After integration  No

Timeout mode  Off
    
```

Select item

Main menu


- Read image
- Find spots
- Edit spots
- Clear spots
- Select images
- Autoindex
- Estimate mosaicity
- Predict

muld1x1_301.mar2000

Min 1 Max 2207 Cursor position

Overlay on Contrast

Colour Black on whi Mag x4



xdl_io_window

Input reply

WARNING, the image display will be closed.
 Do you want to save the current parameters ? (Y)
 Name of save file (muld1x1_301.sav) :refined.sav

This file can be used as input to a processing job.
 Do you want to exit (Y):_

Circles

Beam / backstop

Output

```

Pixel X,Y      2001  0
XC,YC mm      300.0  0.0
Resolution     0.00
Indices        0  0  0
F Phi 0.00 width 0.00
Intensity      0
Sigma          0

Spacing A      0.000
Average        0.0
Rms            0.0
Number         0
Zoomfactor     0
Circle resolution A
      0.0  0.0  0.0  0.0
Phi      279.00 280.00
Missets ThetaX,Y,Z
      0.02  0.10 -0.01
    
```

Blue: fulls, Yellow: partials, Red: overlaps
 Green: too wide in phi

Running the STRATEGY option

- ❑ In the GUI, click on the STRATEGY button
 - ❑ On the command-line, type STRATEGY.
 - ❑ In either case, determining a suitable strategy for data collection once you have a cell, orientation and crystal symmetry is straightforward.
 - ❑ Then run TESTGEN to check for overlaps
-

Input replyUp  Dn

To test the full rotation for this Laue group just type STRATEGY at the prompt.
To test a smaller rotation (eg 60 degrees) in 2 different segments type:
STRATEGY SEGMENTS 2 ROTATE 60
A suitable speedup factor will be set by default but can be overridden
eg STRATEGY ROTATE 60 SEGMENTS 2 SPEEDUP 10
To try specific phi ranges type : STRATEGY START -20 END 60

Do you want to proceed (Y):

MOSFLM =>strategy

Running strategy in default mode. Type GO to continue or ABORT to stop.

MOSFLM =>go

SPEEDUP factor set to 15.0

Angles between a axis and X,Y,Z axes 54.552 35.492 88.463

Angles between b axis and X,Y,Z axes 174.303 84.515 91.537

Angles between c axis and X,Y,Z axes 91.615 90.736 1.775

The c axis is closest to the rotation axis (angle 1.77)

Rotation angle to get the a axis in XZ plane -54.54

Rotation angle to get the b axis in XZ plane 5.49

Unique axis is: c

Rotation angle to get the unique axis into the YZ plane: 245.5

***** WARNING *****

***** WARNING *****

***** WARNING *****

***** WARNING *****

With the crystal in the current orientation it is not possible to collect 100% of the data, some will be lost in the cusp.

To avoid a cusp, the c-axis must be at an angle between 74.5 and 15.5 degrees away from the rotation axis (Z)

You may want to adjust the orientation with the goniometer arcs.

Generating reflection list

Generated reflections for phi range -145.0 to -25.0

Sorting the generated list.

Generating the list of unique reflections.

Merging the generated and unique reflection list.

Checking completeness of data

Optimum rotation gives 97.7% of unique data.

This corresponds to the following rotation range(s):

From 255.0 to 315.0 degrees

Type "STATS" at prompt for full statistics.

It may be possible to get a higher completeness using two segments, or a slightly greater rotation angle.

Use the TESTGEN option to determine appropriate oscillation angles to avoid spatial overlaps.


STRATEGY =>_

Checking the output (1)


- There are two useful log files;
 - SUMMARY - this is of most use when viewed with the CCP4 graph viewer LOGGRAPH, as it contains graphs of parameters which have varied through the data processing.

 - mosflm.lp - this can be very large, and contains a complete record of the experiment.
-


Loggraph summary.1



Loggraph summary.1



Loggraph summary.1



Checking the output (2)

- If everything has gone right so far;
 - check the MTZ file - is it as you'd expect?

More on Processing with Mosflm

- what do you do about the warnings?
 - what if are real problems?
 - the new ccp4i Mosflm task
-

HEADER INFORMATION FROM OUTPUT MTZ FILE

Logical Name : muldlx1_301.mtz Filename: muldlx1_301.mtz

<snip>

* Number of Columns = 18

* Number of Reflections = 43904

* Missing value set to NaN in input mtz file

* Number of Batches = 22

* Column Labels :

H K LM/ISYM BATCH I SIGI IPR SIGIPR FRACTIONCALC XDET YDET...

* Column Types :

H H H Y B J Q J Q R R R R R R I I R

*** For information only. ***

PARTIALS INCLUDED IN POSITIONAL REFINEMENT AND PROFILES

=====

Because there were rather few fully recorded reflections...

<snip>

*** Warning messages ***

TANGENTIAL OFFSET UNSTABLE

=====

The tangential offset parameter (TOFF) is varying more...

<snip>

SPOT OVERLAP

=====

Adjacent spots overlap. This will produce systematic errors...

<snip>

EXCESSIVE NUMBER OF BADSPOTS

=====

At least some images have rather a lot of badspots...

<snip>

TOO MANY BACKGROUND PIXELS OVERLAPPED BY NEIGHBOURING SPOTS

=====

For some of the standard profiles, more than half the backgr...

*** Warning messages ***

TANGENTIAL OFFSET UNSTABLE

=====

The tangential offset parameter (TOFF) is varying more than it should. (Maximum variation is 0.15 mm)


If there are large changes in both TOFF and ROFF or CCOMEGA, this suggests that the refinement is unstable.

In this case, it is best to fix the TOFF parameter:

REFINEMENT FIX TOFF

If known the correct value can be input:eg

DISTORTION TOFF 0.17 If not known, the mean refined value can be used. In such cases ROFF should also be FIXED.



Tables in File

Refined detector parameters

Post refinement

Graphs in Selected Table

Camera constants CCX CCY (mm) CCOM (deg) v image

Distance (mm) v image

Yscale v image

Tilt and twist (hundredths of deg) v image

ROFF and TOFF (mm) v image

SPOT OVERLAP

=====

<snip>

The minimum allowed spot separation (SEPARATION keyword) was 1.6 1.6mm. The actual spot size determined by the mask optimisation is 1.7 by 1.4mm in the centre of the image and the largest spot size is 2.8 by 2.8mm.

The separation given should be at least as large as the spot size in the centre of the image (keyword SEPARATION).

Check standard profiles carefully to ensure that the optimisation of the raster parameters has worked correctly.

The effective size of the spots can be controlled by PROFILE TOLERANCE keywords. If the peak regions look too large (ie they include too much of the tails of the spot), try increasing TOLERANCE (current value 0.010) by eg 0.005 and see if profiles look better. (Increasing TOLERANCE will decrease spot size).

As a last resort the profile optimisation can be turned off using keywords PROFILE NOOPT.

In cases of serious overlap, (ie if the pattern is very dense), then the SEPARATION CLOSE option should be used (eg SEPARATION 1.0 1.0 CLOSE) and it may also help to suppress profile optimisation in these cases (PROFILE NOOPT) keyword. See help library for details.

**** WARNING **** Peaks of neighbouring reflections overlap, data quality will be impaired

<snip>

Profile for box 1

X limits 0 to 92 mm, Y limits 0 to 92 mm

Number of reflections in profile 299 RMSBG 9.5 Profile factor 0.50

```
-0-]-0-0*0*]*0*0*0*0*]*0*]*0*0*0*0*]*]*]*]*]*]*]-]
-0-0-0-0*0*0*0*0*0*0*0*0*0*0*1*0*]*0*0*]*]*]*]*]*]*]-]
-0-]-]-0*1*1*1*1*1*0*0*0*0*1*1*0*0*0*0*0*]*]*]*]*]-]-]
-0-0-0-0-0*0*0*0*0*1*2*2*1*1*2*1*0*0*1*0*]*]*]*]-]-]-]
-1-0-0-0-0-0*1*1*1*2*2*1*1*2*1*1*1*1*0*1*0-]-]-]-]
-0-]-0-1-1-1-1 2 2 1 2 3 3 4 3 3 2 2*1*1-0-0-0-]-]
-0-0-0-1-2-1 1 2 3 3 4 5 6 7 6 5 4 2 1-0-]-]-0-0-]
-1-1-1-1-1 1 2 3 4 5 7 9 BC B 8 6 4 2 1-1-0-0-0-0
-2-1-1-1 1 2 2 4 5 8 B F I I G C 9 6 4 2 1-0-0-]-0
-1-1-1 1 1 2 3 4 7 B H M Q P M H B 8 5 3 2 0-0-0-0
-1-2-1 2 2 2 2 5 9 G N U X W R L F 9 5 2 1 0-0-0-0
-0-1-1 2 2 3 3 6 C J Q W Z X S L E 9 5 2 1 1-1-0-0
-1-2-2 3 2 2 4 6 C J P V W U Q J D 8 5 4 1 1-0-0-0
-1-2-2 2 2 2 4 6 A G L P Q O K F A 7 5 3 2 1-0-1-1
-1-1-1-3 2 3 3 5 8 C F H H G D A 7 5 3 2 1-1-0-1-1
-1-1-1-2-2 2 3 4 6 8 A B A 9 8 6 5 3 3 1-0-1-1-0-1
-0-0-0-1-1-2 2 3 4 5 6 6 7 6 5 4 4 3 3-2-1-1-2-]-0
-]-1-1-1-2*2*2 3 3 4 4 4 5 4 4 3 3 2-2-1-1-1-1-0-0
-]-]-1-2*2*2*2*2*3*4*3*3*3*3*3*2*2*1*2-2-1-1-1-1-0
-]-]-0*1*2*2*1*2*2*2*3*2*2*2*2*2*3*3*1*2-1-]-0-0-1
-]-]*]*1*2*2*2*2*1*2*2*3*2*2*1*2*2*2*1*1*0-1-0-0-1
-]*]*]*]*2*2*2*1*1*2*2*2*2*2*1*1*1*1*2*0*0-1-0-0-0
-]*]*]*]*]*1*1*2*1*1*1*2*2*2*1*1*2*2*1*0*0-1-0-]-0
```

EXCESSIVE NUMBER OF BADSPOTS

=====

At least some images have rather a lot of badspots (Maximum number 42). They are rejected on the basis of:

- 1) Poor profile fit (PKRATIO >3, controlled by REJECTION PKRATIO). 6
- 2) Too large a BGRATIO (too much background variation, controlled by REJECTION BGRATIO).
- 3) Too large a background gradient (controlled by REJECTION GRADMAX) 578
- 4) Intensity negative and more than 5 sigma. 38

Look at the list of badspots to see what category they fall under.

Poor profile fit is often the result of changes in ROFF, TOFF or CCOMEGA between successive images when using the ADDPART option.

Very intense images can have unusually large gradients, GRADMAX may have to be changed from the default


A pixel dump of the BADSPOTS can be obtained using REJECTION PLOT if the reason for their rejection is not clear

TOO MANY BACKGROUND PIXELS OVERLAPPED BY NEIGHBOURING SPOTS

=====
For some of the standard profiles, more than half the background pixels are flagged as being overlapped by neighbouring spots (in the worst case, 52.2 % are overlapped). You should use the SEPARATION CLOSE keywords, eg
SEPARATION 1.0 1.0 CLOSE

More on Processing with Mosflm

- what do you do about the warnings?
 - what if there are real problems?
 - the new ccp4i Mosflm task
-



12	103	mC	43.79	219.23	43.93	90.0	90.2	79.6	C2
11	103	oC	43.79	219.23	43.93	90.0	90.2	100.4	C222, C2221
10	7	tP	43.79	43.93	107.84	90.0	91.1	90.2	P4, P41, P42, P43, P422, P4212
9	6	mC	61.93	62.13	107.84	89.2	90.8	89.8	C2
8	6	mC	62.13	61.93	107.84	90.8	90.8	90.2	C2
7	6	oP	43.79	43.93	107.84	90.0	91.1	90.2	P222, P2221, P21212, P212121
6	6	mP	43.79	107.84	43.93	90.0	90.2	91.1	P2, P21
5	6	oC	61.93	62.13	107.84	89.2	90.8	89.8	C222, C2221
4	1	mP	43.79	43.93	107.84	90.0	91.1	90.2	P2, P21
3	0	mP	43.79	43.93	107.84	90.0	91.1	90.2	P2, P21
2	0	aP	43.79	43.93	107.84	90.0	88.9	89.8	P1
1	0	aP	43.79	43.93	107.84	90.0	91.1	90.2	P1

Select a solution AND a spacegroup from list above (eg 3 p42) or 0 to abandon or T to cha

The solution and direct beam position will now be refined; reflections which deviate by m
 than the sigma cutoff from their calculated position will be excluded from the refinement

Positional sigma cutoff [2.50]:

Refining solution #10 with P4 (number 75) symmetry imposed

Using 322 indexed reflections (out of 366 spots found, {delta(XY) <= 2.5 sigma}),
 final sd in spot positions is 0.38mm and in phi 0.78 degrees

Refined cell parameters 44.03 44.03 107.87 90.00 90.00 90.00


Do you want to update cell parameters (Y):

Beam coordinates of 150.30 150.00 have been refined to 150.82 149.63

***** WARNING ***** WARNING ***** WARNING ***** WARNING ***** WARNING

This is a shift of 0.64mm or 0.298 times the minimum spot separation of ca 2.14mm.
 Do you want to accept the new direct beam position? (answer Y or N!) :y

Do you want to accept this solution (Y) :_



Cell refinement is complete

Starting cell	44.033	44.033	107.873	90.000	90.000	90.000
Refined cell	43.953	43.953	107.685	90.000	90.000	90.000

Rms positional error (mm) as a function of cycle for each image.

Image	1	2	3	88	89	90
Cycle 1	0.160	0.168	0.181	0.738	0.663	2.608

YSCALE as a function of cycle for each image:

Image	1	2	3	88	89	90
Cycle 1	1.000	1.002	1.000	0.966	0.964	0.970

Detector distance as a function of cycle for each image:

Image	1	2	3	88	89	90
Cycle 1	149.2	149.1	150.1	151.3	151.0	151.6

Refined mosaic spread (excluding safety factor): 0.12

Missets for first image (1) -0.03 0.02 -0.02

Missets for last image (90) 0.06 -0.27 -0.03

The current missets are for the last image to be processed.

If you want to integrate the data starting at the first image, you should reset the missetting angles.

Reset missets to those of the first image ? (Y)_

12	103	mC	43.79	219.23	43.93	90.0	90.2	79.6	C2
11	103	oC	43.79	219.23	43.93	90.0	90.2	100.4	C222, C2221
10	7	tP	43.79	43.93	107.84	90.0	91.1	90.2	P4, P41, P42, P43, P422, P4212
9	6	mC	61.93	62.13	107.84	89.2	90.8	89.8	C2
8	6	mC	62.13	61.93	107.84	90.8	90.8	90.2	C2
7	6	oP	43.79	43.93	107.84	90.0	91.1	90.2	P222, P2221, P21212, P212121
6	6	mP	43.79	107.84	43.93	90.0	90.2	91.1	P2, P21
5	6	oC	61.93	62.13	107.84	89.2	90.8	89.8	C222, C2221
4	1	mP	43.79	43.93	107.84	90.0	91.1	90.2	P2, P21
3	0	mP	43.79	43.93	107.84	90.0	91.1	90.2	P2, P21
2	0	aP	43.79	43.93	107.84	90.0	88.9	89.8	P1
1	0	aP	43.79	43.93	107.84	90.0	91.1	90.2	P1

Select a solution AND a spacegroup from list above (eg 3 p42) or 0 to abandon or T to change

The solution and direct beam position will now be refined; reflections which deviate by more than the sigma cutoff from their calculated position will be excluded from the refinement

Positional sigma cutoff [2.50]:

Refining solution # 7 with P222 (number 16) symmetry imposed

Using 322 indexed reflections (out of 366 spots found, {delta(XY) <= 2.5 sigma}), final sd in spot positions is 0.22mm and in phi 0.64 degrees

Refined cell parameters 45.77 43.87 107.79 90.00 90.00 90.00

Do you want to update cell parameters (Y):

Beam coordinates of 150.30 150.00 have been refined to 150.61 149.60

This is a shift of 0.51mm or 0.100 times the minimum spot separation of ca 5.05mm. Do you want to accept the new direct beam position? (Y) :

Do you want to accept this solution (Y) :_

Cell refinement is complete

Starting cell	43.864	45.758	107.796	90.000	90.000	90.000
Refined cell	43.962	45.897	108.009	90.000	90.000	90.000

Rms positional error (mm) as a function of cycle for each image.

Image	1	2	3	88	89	90
Cycle 1	0.066	0.062	0.061	0.256	0.305	0.190
Cycle 2	0.112	0.148	0.110	0.231	0.274	0.220

YSCALE as a function of cycle for each image:

Image	1	2	3	88	89	90
Cycle 1	1.000	1.000	1.000	1.000	0.996	0.997
Cycle 2	1.000	0.999	1.000	0.997	0.999	0.999

Detector distance as a function of cycle for each image:

Image	1	2	3	88	89	90
Cycle 1	150.0	150.0	150.1	149.8	150.4	150.1
Cycle 2	150.5	150.6	150.5	150.9	150.6	150.6

Refined mosaic spread (excluding safety factor): 1.60

Missets for first image (1)	0.06	-0.01	0.15
Missets for last image (90)	0.17	-0.33	0.21

The current missets are for the last image to be processed.

If you want to integrate the data starting at the first image, you should reset the missetting angles.

Reset missets to those of the first image ? (Y)_

More on Processing with Mosflm

- what do you do about the warnings?
 - what if there are real problems?
 - the new ccp4i Mosflm task
-

Enter job title (use only alphanumeric, spaces, brackets and underscores)

Help

*Run MOSFLM in batch mode to integrate diffraction images*Job title

Take parameters from command file

Specify where the images are located Use current working directory /Users/harry/test/muldlx1Image file name template: Matrix file /Users/harry/test/muldlx1/muldlx1_301.mat

Browse

View

First image Last image MTZ out

Browse

View

Dataset and Harvesting Parameters Create harvest file in project harvesting directoryHarvest project name and dataset name Crystal name for this dataset *Crystal Parameters*Space group Cell a b c alpha beta gamma Mosaicity *Detector Parameters* Set detector type to and over-ride default determined by MOSFLMDetector to crystal distance mm Detector swing angle degreesMinimum spot separation: x y mm Close spots? Detector gain Scanner ADC offset Overload cutoff value Overload flagging value Measurement box: Width Height Corner cut-off Horiz. bg Vert. bg *Distortion Parameters**Main Beam Position*Beam position: x y (mm) Beam swung out

Run

Save or Restore

Close

Enter job title (use only alphanumeric, spaces, brackets and underscores)

Help

*Run MOSFLM in batch mode to integrate diffraction images*Job title

Take parameters from command file

Specify where the images are located Use current working directory /Users/harry/test/muldlx1Image file name template: Matrix file /Users/harry/test/muldlx1/muldlx1_301.mat

Browse

View

First image Last image MTZ out muldlx1_0.mtz

Browse

View

Dataset and Harvesting Parameters Create harvest file in project harvesting directoryHarvest project name and dataset name Crystal name for this dataset *Crystal Parameters*Space group

Cell a 125.7409 b 125.7409 c 74.0908 alpha 90.0000 beta 90.0000 gamma 120.0000

Mosaicity *Detector Parameters* Set detector type to and over-ride default determined by MOSFLMDetector to crystal distance mm Detector swing angle degreesMinimum spot separation: x y mm Close spots? Detector gain Scanner ADC offset Overload cutoff value Overload flagging value Measurement box: Width Height Corner cut-off Horiz. bg Vert. bg *Distortion Parameters**Main Beam Position*Beam position: x y (mm) Beam swung out

Run

Save or Restore

Close

Enter job title (use only alphanumeric, spaces, brackets and underscores)

Help

*Run MOSFLM in batch mode to integrate diffraction images*Job title

Take parameters from command file

Specify where the images are located Use current working directory /Users/harry/test/muldlx1Image file name template: Matrix file /Users/harry/test/muldlx1/muldlx1_301.mat

Browse

View

First image Last image MTZ out

Browse

View

Dataset and Harvesting Parameters Create harvest file in project harvesting directoryHarvest project name and dataset name Crystal name for this dataset *Crystal Parameters*Space group Cell a b c alpha beta gamma Mosaicity *Detector Parameters* Set detector type to and over-ride default determined by MOSFLMDetector to crystal distance mm Detector swing angle degreesMinimum spot separation: x y mm Close spots? Detector gain Scanner ADC offset Overload cutoff value Overload flagging value Measurement box: Width Height Corner cut-off Horiz. bg Vert. bg *Distortion Parameters**Main Beam Position*Beam position: x y (mm) Beam swung out

Run

Save or Restore

Close

Data Collection & Data Analysis

Scaling and merging with Scala

Gwyndaf Evans

Diamond Light Source, UK.

Contents

- ❑ sources of systematic errors
 - ❑ parameterization of scaling
 - ❑ estimates of data quality
 - ❑ estimation of standard errors
 - ❑ outlier rejection
-

Introduction

- scaling and merging is the most important diagnostic step in terms of data quality.
 - it is important that it be performed as soon as possible after data collection and preferably during.
 - in many cases it is a straightforward procedure but can become complicated. An understanding of the underlying principles is important.
-

Steps in scaling

- Choose scaling model
 - Should reflect the experiment
 - X-ray source intensity variation
 - Changes in diffracting volume
 - Sample or air absorption
 - Radiation damage
 - Analyze results
 - Should the sample be discarded?
 - What is the useful resolution?
 - Are there outliers or bad images?
 - What is the spacegroup?
-

Factors affecting scale

Incident beam

- intensity
- size
- primary beam absorption

Sample

- diffracting volume
- diffracted beam absorption

Detector

- calibration
- time stability
- bad pixels

Miscellaneous

- beam stop and cryo-stream shadows
-

Incident beam related factors

□ Synchrotron

- smooth decay of beam intensity
 - any discontinuities (e.g. beam injection) should be noted and included in scaling model
 - illuminated volume
 - shutter synchronization/goniometer rotation speed
-

Crystal related factors

□ Sample absorption

- diffracted beam absorption (shape dependent)
- important for weak anomalous signal

□ Radiation damage

- can be significant on high brilliance sources
 - difficult to correct for
 - modeled as change in relative B-factor
 - extrapolation to zero dose
-

Detector related factors

□ calibration errors

- spatial distortion
 - non-uniformity of response
 - time stability
 - bad pixels
-

Miscellaneous factors

unavoidable

- zingers

avoidable

- beam stop shadow
 - cryo-stream shadow
 - should be dealt with at integration stage
-

Determination of scale factors

What information do we have?

Scales are determined by comparison of symmetry-related reflections, i.e. by adjusting scale factors to get the best internal consistency of intensities. Note that we do not know the true intensities and an internally-consistent dataset is not necessarily correct. Systematic errors will remain

$$\text{Minimize } \Delta\Phi = \sum_{hl} w_{hl} (I_{hl} - 1/k_{hl} \langle I_h \rangle)^2$$

I_{hl} l 'th intensity observation of reflection \mathbf{h}

k_{hl} scale factor for I_{hl} $\langle I_h \rangle$ current estimate of I_h

$g_{hl} = 1/k_{hl}$ is a function of the parameters of the scaling model

$$g_{hl} = g(\square \text{ rotation/image number}) \cdot g(\text{time}) \cdot g(\mathbf{s})$$

...other factors

Primary beam \mathbf{s}_0

B-factor Absorption

Parametrization



$$g_{hl} = g(\text{rotation/image number}) \cdot g(\text{time}) \cdot g(\mathbf{s}) \cdot \dots \text{other factors}$$

Primary beam s_0

B-factor

Absorption

eg "tails"





scale is smooth function of spindle rotation *or* discontinuous function of image (batch) number (usually less appropriate)

$$g(\text{time}) = \exp[+2B(\text{time}) \sin^2\theta/\lambda^2]$$


essentially a time-dependent radiation damage correction

Absorption correction

- ❑ serious for longer wavelengths and larger irregular samples
- ❑ we typically do not have enough data to find true absorption corrections
- ❑ good corrections require data recorded with sample in multiple orientations
- ❑ some strategies can minimize effect of absorption on quality of observed anomalous differences
- ❑ inverse beam
- ❑ simultaneous measurement of I^+ & I^-




Spherical harmonics





$$A(\mathbf{s}) = \sum_{lm} c_{lm} Y_{lm}(\theta_p, \phi_p)$$

- linear coefficients c_{lm} determined as parameters.
- Note the surface is not centrosymmetric (see e.g. equator $\theta = 90^\circ$)
- i.e. different corrections are applied to I+ & I-
- Graph along lines of latitude
- Dashed lines where there is no data






Absorption: sample dataset

- ❑ Rotating anode (RU200, Osmic mirrors, Mar345)
- ❑ 100 images, 1° , 5 min/ $^\circ$, resolution 1.8Å
- ❑ Scaled with and without secondary beam correction (sub-keyword: absorption)
- ❑ Secondary beam correction (absorption) improves the data



Absorption correction and phasing



□ Phasing (SHARP) is better with absorption correction even after solvent flattening

SCALA 'scales' options

- Typical scaling options are

scales rotation spacing 5 absorption 6 bfactor on brotation spacing 20

tie surface 0.005

restrain secondary beam surface

maybe relax (eg 0.05) for high absorption

link surface all

same surface for all runs

tie bfactor 0.3

stabilise B-factors

- When to use other options (more complicated cases):-
 - data collected from more than one crystal
 - use different absorption surfaces
 - very low resolution - don't use B-factor
-

How well are the scales determined?

- Note that determination of scaling parameters depends on symmetry-related observations having **different** scales. If all observations of a reflection have the same value of the scale component, then there is no information about that component and it remains as a systematic error in the merged data (this may well be the case for absorption for instance)
 - Thus to get intensities with the lowest absolute error, the symmetry-related observations should be measured in as different way as possible (eg rotation about multiple axes). This will increase R_{merge} , but improve the estimate of $\langle I \rangle$.
 - Conversely, to measure the most accurate differences for phasing (anomalous or dispersive), observations should be measured in as similar way as possible
-

Scaling datasets together

- For multiple-wavelength datasets, it is best to scale all wavelengths together simultaneously. This is then a *local* scaling to minimize the difference between datasets, reducing the systematic error in the anomalous and dispersive differences which are used for phasing
 - Other advantages of simultaneous scaling:-
 - rejection of outliers with much higher reliability because of higher multiplicity
 - correlations between ΔF_{anom} and ΔF_{disp} indicate the reliability of the phasing signal
 - approximate determination of relative f'' and relative f' values
 - *In Scala, this is now automatic if multiple datasets are present in the input file*
-

Results of scaling: R-factors


- (a) $R_{\text{merge}} (R_{\text{sym}}) = \frac{\sum | I_h | - \langle I_h \rangle}{\sum \langle I_h \rangle}$
 - This is the traditional measure of agreement, but it increases with higher multiplicity even though the merged data is better
 - (b) $R_{\text{meas}} = R_{\text{r.i.m.}} = \frac{\sum (n/n-1) | I_h | - \langle I_h \rangle}{\sum \langle I_h \rangle}$
 - The multiplicity-weight R-factor allows for the improvement in data with higher multiplicity. This is particularly useful when comparing different possible point-groups
 - Diederichs & Karplus, *NSB.*, 4, 269-275 (1997)
-

Results of scaling: Intensities and standard deviations

□ Scala compares the estimated standard deviation $\sigma(I)$ to the observed scatter, and tries to correct $\sigma(I)$ by a multiplication factor. This is done using a *normal probability plot*. A correction as a function of intensity is also done, but this is not yet automatic

□ $\sigma(I)' = \text{Sdfac} * \text{Sqrt} [\sigma^2(I) + (\text{Sdadd} * I)^2]$

□ The corrected $\sigma(I)$ is compared with the intensities: the most useful statistic is $\langle I \rangle / \sigma(\langle I \rangle)$ (labelled Mn(I)/sd in table)




$\langle I \rangle / \sigma(\langle I \rangle) > .gt. \sim 2$

Estimation of errors



- keyword for estimation of standard errors
 - `sdcorrection sdfac 1.5 sdadd 0.015`

Normal probability plots



Results of scaling: Correlation coefficients

- CC's between different estimates of anomalous or dispersive differences (ΔI) indicate reliability of the phasing signal




Correlations between half-sets
(data randomly halved)

Correlations between $(I^+ - I^-)$
for MAD data


Results of scaling: Common sense checks

- Are some parts of the data bad?
- Analysis of Rmerge against batch number gives a very clear indication of problems local to some regions of the data. Perhaps something has gone wrong with the integration step, or there are some bad images




- Here the beginning of the dataset is wrong due to problems in integration (e.g. poor orientation matrix in MOSFLM at start of job.)

Results of scaling: Physical sense checks





- These scale factors follow a reasonable absorption curve



- These B-factors are not sensible (fluctuations may be damped by "tie bfactor")

Spacegroup determination



- Lattice type: hexagonal (from indexing)
- Merge in different Laue groups
- Point-group P321 is highest symmetry with good merging (Rmeas)

- Systematic absences
 - $l = 3n$
- Spacegroup = P3121
 - could be P3221

Outliers

- Some (a small proportion) of observations may be rejected as being improbable
 - Unreasonably large
 - default limit in Scala $E > 10$
 - removes **some** ice|salt spots & zingers
 - Unreasonable disagreement with other observations
 - easier with high multiplicity
 - not clear what to do between *e.g.* I+ & I- where there is real difference
-

Outliers: possible causes

- outside reliable area of detector (*e.g.* behind shadow)
 - specify backstop shadow, calibrate detector
 - ice spots
 - do not get ice on your crystal!
 - zingers
 - bad prediction (spot not there)
 - improve prediction
 - spot overlap
 - lower mosaicity, smaller slice, move detector back
 - deconvolute overlaps
 - multiple lattices
 - find single crystal
-

Outliers: detection

- ❑ Detection of outliers is easiest if the multiplicity is high
- ❑ Removal of spots behind the backstop shadow does not work well at present: usually it rejects all the good ones, so **tell Mosflm where the backstop shadow is**
- ❑ Inspect the ROGUES file to see what is being rejected (at least occasionally)

The ROGUES file contains all rejected reflections (flag "*", "@" for I+- rejects, "#" for Emax rejects)

TotFrc = total fraction, fulls (f) or partials (p)

Flag I+ or I- for Bijvoet classes

$\text{DelI/sd} = (\text{Ihl} - \text{Mn(I)others}) / \sqrt{[\text{sd}(\text{Ihl})^2 + \text{sd}(\text{Mn(I)})^2]}$

h	k	l	h	k	l	Batch	I	sigI	E	TotFrc	Flag	Scale	LP	DelI/sd	d(A)
Xdet	Ydet	Phi													
(measured)	(unique)														
-2	-2	0	2	2	0	1220	24941	2756	1.03	0.95p	I-	2.434	0.031	-1.1	30.40
1263.7	1103.2	210.8													
-4	2	0	2	2	0	1146	9400	2101	0.63	0.99p	*I+	3.017	0.032	-6.7	30.40
1266.4	1123.3	151.3													
4	-2	0	2	2	0	1148	27521	2972	1.08	1.09p	I-	2.882	0.032	0.0	30.40
1058.8	1130.0	153.2													
2	-4	0	2	2	0	1075	29967	2865	1.13	0.92p	I+	2.706	0.032	1.1	30.40
1060.9	1106.6	94.4													

Acknowledgments

mosflm Andrew Leslie and Harry Powel

scala Phil Evans